


# RevoLix<sup>jr</sup>.

Your Table top surgical Thulium laser

## Applications

ENT  
Pneumology  
Neurosurgery\*  
Gynaecology  
Gastroenterology  
Urology  
Visceral Surgery  
Laparoscopy



### RevoLix jr. – powerful thulium laser in compact design

This compact table top surgical laser is available for soft tissue surgery and unifies all advantageous properties of existing laser principles in a single easy transportable unit:

**RevoLix jr.** offers the cutting and ablation properties like the well known CO<sub>2</sub> laser from a flexible fibre - no need for an articulated mirror arm or hollow waveguide.

**RevoLix jr.** laser radiation achieves excellent haemostasis like Nd:YAG and diode lasers – without the risk of deep tissue penetration.

**RevoLix jr.** tissue effect is shallow like the Holmium laser – however there is no trauma. Cutting edges are smooth and clean. In open surgery there is no splattering.

**RevoLix jr.** laser radiation is delivered to the surgical site by flexible fibres – ideal for endoscopic, laparoscopic and minimally invasive surgery.

### RevoLix jr. – why Two Micron?


The **RevoLix jr.** wavelength at 2.0 micron is almost identical to the well established Holmium laser. This wavelength is known for its suitability for resection and ablation, safe application in an aqueous medium and it generates excellent haemostasis. These excellent properties are due to the efficient absorption of 2.0 micron laser radiation by water which is ubiquitous in any tissue. At the **RevoLix jr.** wavelength of 2.0 micron the absorption is 2.5 times stronger than at the Holmium wavelength providing even more precise cutting in soft tissue surgery.

### RevoLix jr. – why Two Micron cw laser?

In soft tissue surgery efficient and even cutting combined with strong haemostasis is required. This is provided by the **RevoLix jr.** laser in an unmatched manner. The new **RevoLix jr.** laser operates in a continuous wave mode and cuts and vaporizes circulated and white tissue without deep penetration or uncontrolled necrosis. Clean cuts and excellent haemostasis are achieved by moving the beam across the surgical site. Generous Two Micron laser power capacity allows high vaporization and resection rates at no bleeding and short theatre time.

### What are the advantages?

Cutting efficiency of soft tissue and haemostasis is superior to any known alternative. Tissue effect is restricted to less than 1 millimetre beneath the cut. **RevoLix jr.** preserves excellent vision to the surgical site. There is no vision impairment due to bleeding, excessive bubble formation, tissue fragments or ruptured tissue. Unlike visible lasers (KTP, LBO and some diode wavelengths) the surgical site is free of visible glare. The colour neutral laser safety eyewear does not generate discolouration like KTP. Endoscope lenses remain free from splatter when used in open surgery or in laparoscopy.




### RevoLix jr. – why is it safe?

In an aqueous medium the laser effect to tissue is restricted to less than 2 millimetres in front of the tip of the fibre. Any tissue further away is shielded by the medium. The same mechanism is protecting tissue and organs adjacent to the cut. Any tissue more distant than 2 millimeters beneath the cut is unaffected by the **RevoLix jr.** laser. Opposed to KTP, LBO, Nd:YAG and diode lasers this property eliminates the risk of unintentional tissue damage during laser surgery.

### What are the applications?

The **RevoLix jr.** laser system has demonstrated its superior performance in surgical disciplines such as ENT, Pneumology, Neurosurgery\*, Gynaecology, Gastroenterology, Urology, Visceral Surgery and Laparoscopy.



## RevoLix jr. – delivery systems

LISA offers a wide range of specialised delivery systems. Please refer to the **Surgical Laser Accessories** and to the **Surgical Laser Fibres** brochure.

All fibres match with the surgical instrumentation recommended and approved by LISA laser products.

## How does RevoLix jr. match with your theatre ?

The new laser is extremely theatre friendly **RevoLix jr.** does not cause noise strain to theatre personnel. The table top laser fits into any video tower. A laser trolley is available on option. Regular theatre utilities are sufficient (standard power outlet 100 - 240 VAC, no cooling water required). The laser has proven its extreme sturdiness during routine transportation between theatres and other departments.

## RevoLix jr. – operation of the laser

The new laser is easy and safe to operate. The display is easy to understand. Operational modes such as continuous or chopped are selected by push buttons. The power setting is adjusted with a large adjustment wheel. Large characters display the settings on the easy to read display.

The laser is activated by the low profile **Kix** footswitch.



**Kix** footswitch  
Order No. 101 660 043

## Urology


Vaporisation and excision of tumours  
UUT-TCC, UTUC  
Bladder neck incision  
Opening of strictures  
Partial nephrectomy  
Laparoscopy



## Haemorrhoids

## Gynaecology

Excision of polyps  
Endometriosis  
Hysterectomy  
Adhesiolysis  
Conisation  
Condylomata  
Myomectomy




## Neurosurgery\*

Fenestration of cysts  
Ventriculocysternosomy  
Catheter recovery  
3rd ventriculostomy  
Tumour resection  
Haemostasis



## ENT

Excision of tumours  
Excision of granulomas  
Tonsillectomy  
UVPP



## Pneumology

Bronchoscopy  
Airway recanalization  
Desobstruction  
Tissue coagulation



## General Surgery



## Technical Specifications

RevoLix<sup>jr</sup>.



RevoLix<sup>jr</sup>. 30

<b>Laser system</b>	continuous wave DPSS laser	continuous wave DPSS laser
<b>Wavelength</b>	2013 nm	2013 nm
<b>Power at fibre tip</b>	1 - 15 Watt cw (adjustable)	1 - 30 Watt cw (adjustable)
<b>Chopped mode</b>	50 - 1000 ms	50 - 1000 ms
<b>Frequency</b>	0.5 - 10 Hz	0.5 - 10 Hz
<b>Beam delivery</b>	wide range of fibres	wide range of fibres
<b>Aiming beam</b>	635 nm (red) or 532 nm, (green) max. 1.3 mW (adjustable), laser class 3R	
<b>Mains supply</b>	100 - 240 VAC 50/60 Hz, 6 A max.	100 - 240 VAC 50/60 Hz, 6 A max.
<b>Cooling system</b>	Air cooling	Air cooling
<b>Dimensions</b>	H 160 x W 460 x L 630 mm	H 160 x W 460 x L 630 mm
<b>Weight</b>	20 kg	20 kg
<b>Environmental conditions</b>	15 - 28 °C / 10 - 90 % humidity (non-condensing)	

Safety Standards: IEC 60601

CE acc. Council Directive 93/42/EEC

U.S. federal law restricts this device to sale by or on the order of a physician.

**Laser Trolley**

Optional laser trolley



### IMPORTANT NOTICE:

The information provided is a general overview of potential clinical applications of the described products. National health care regulations vary between countries and may exclude certain clinical applications at your location. The user assumes responsibility to be updated about national deviations from the applications listed above.

\* In the USA the products are not intended for use in clinical applications in neurosurgery.

U.S. federal law restricts these devices to sale by or on the order of a physician

Specifications are subject to change without notice.

Made in Germany 2019 - 02

Brochure Order no. 036 004 110

VISIBLE AND INVISIBLE LASER RADIATION  
AVOID EYE OR SKIN EXPOSURE TO  
DIRECT OR SCATTERED RADIATION  
CLASS 4 LASER PRODUCT

IEC 60825-1:2007



DPSS LASER:  $\lambda = 2013 \text{ nm}$  /  $P_e = 36 \text{ W} / \text{cw}$

DIODE LASER:  $\lambda = 532 \text{ nm}$  /  $P_e = 5 \text{ mW} / \text{cw}$



LISA Laser Products GmbH  
Albert-Einstein-Str. 4  
37191 Katlenburg-Lindau  
Germany  
fon: +49 5556 9938-0  
fax: +49 5556 9938-10  
[info@lisalaser.de](mailto:info@lisalaser.de) [www.lisalaser.de](http://www.lisalaser.de)

**LISALaser**  
PRODUCTS  
An OmniGuide Company