



The Way Aesthetic Treatment Should Be





THE WAY AESTHETIC TREATMENT SHOULD BE

JetPro is the latest JetPeel® device, using jet stream technology for fast, non-invasive transdermal delivery and visibly fresh rejuvenated skin. The new JetPro device is quiet and sleek, delivering the world's most beautifully comfortable and effective treatment.

JetPro harnesses advanced hydra dermabrasion technology to provide deep transdermal treatment with no needles, no skin breakage, no scar tissue and no recovery period. Just excellent results and a pain-free, relaxing treatment so all clients can reveal the very best version of themselves, inside and out.





SPECIFICATIONS

Dimensions:
460mm x 560mm x 680mm

Weight:
27 kg

Mobile:
Integral base-wheels



Max. Source Pressure:
90 psi

Liquid Source:
Disposable solution container

JETPEEL TECHNOLOGY: PERFECTLY COMFORTABLE AND EFFECTIVE

JetPro is based on advanced, exclusive jet propulsion technology adapted by TavTech and developed for the world of aesthetics. It's a revolutionary treatment method that breaks all the barriers of standard beauty care.

TREATMENT SEQUENCE



STEP
01

LYMPHATIC DRAINAGE

A pleasant and relaxing massage useful for detoxifying the skin and stimulating micro-circulation

STEP
02

HYDRA DERMABRASION

Removal of dead skin cells prepares the skin to receive nourishment in the Infusion stage

STEP
03

HYDROPORATION

This step provides effective and deep delivery of active ingredients into the skin

STEP
04

BOOSTER INFUSION

This optional step can be included to boost the treatment result





JETPEEL HANDPIECE THE MIRACLE WAND

The patented JetPeel handpiece is the core state-of-the-art technology that enables delivery of fine misty liquid jet streams at subsonic speeds. The jet streams stretch the skin creating tiny channels that open the way to non-invasive transdermal delivery of nutritive treatment formulas in the deepest skin layers. The JetPeel handpieces facilitate direct and instant delivery for beautiful results you can see already as soon as the treatment is over.

JetPro is compatible with a range of handpieces, designed for various uses and treatments:



MAGIC - Single nozzle, ideal for mesotherapy and anti-wrinkle treatments



SINGLE - Single larger nozzle, for solutions with higher molecular weight



TRIPLE - Three nozzles, for use with most Jet Solutions and treatments

JET SOLUTIONS FOR EVERY AESTHETIC NEED

Jet Solutions are nutritive power treatments that deeply penetrate the skin with the non-invasive, pain-free and comfortable JetPro delivery system. Used in conjunction with JetPro, Jet Solutions offer immediate, visible results for a wide range of aesthetic challenges, including:




ANTI-AGING




RENEWAL



WHITENING



SWEAT



HAIR CARE



Look Beautiful. Feel Amazing.

For professional aestheticians and skincare specialists, JetPro combines the best of both worlds: instant, clinically proven results that look beautiful, in a treatment method that is pain-free and comfortable, leaving the client feeling amazingly fresh and rejuvenated.



No pain

JetPro is completely pain free, with no needles or skin breakage. The non-invasive jet streams penetrate deep into the skin in a totally comfortable, relaxing and enjoyable experience.



No downtime

100% non-invasive transdermal delivery means there is no scar tissue and no recovery time. The client leaves the treatment in full comfort and can immediately return to their daily routine.



Fast visible results

Instant delivery deep into the skin with the JetPeel handpiece means that results are immediate and can be seen already straight away.



Confident & beautiful

JetPro liberates clients from standard aesthetic care, providing a comfortable and relaxing treatment that is instantly effective and rejuvenating - empowering them to be the very best version of themselves.



TAVTECH LTD
T: +972.3.536.0515
E: info@jetpeel.com

40 Haatzmaut St.
Yehud 5630417
Israel



www.jetpeel.com