



Dantec™ Clavis™

Handheld device for EMG/STIM-guided injections

Simple sophistication

Dantec Clavis is designed to be simple and intuitive to provide accurate guidance in minimal time

A fast and easy pathway to injection of medicinal products

The Clavis device is designed to provide the functionality of EMG recording and stimulation, with the mobility and simplicity Clavis can deliver. Clavis is the first device to combine the quality of larger systems with the ease and flexibility of a handheld "smartphone-size" device.

- EMG and STIM functions in one compact device
- Easy to use & operate
- Battery powered
- Perfect match with Bo-ject needles

Disposable Bo-ject® Needle Electrodes

Order #	Size	Lead Length	Instrument End
9013S0422	25 x 0.30 mm	90 cm	1.5 mm female TP
9013S0432	25 x 0.41 mm	90 cm	1.5 mm female TP
9013S0442	37 x 0.41 mm	90 cm	1.5 mm female TP
9013S0472	37 x 0.46 mm	90 cm	1.5 mm female TP
9013S0452	50 x 0.46 mm	90 cm	1.5 mm female TP
9013S0462	75 x 0.64 mm	90 cm	1.5 mm female TP

Each package contains 12 electrodes in individual pouches

Pre-gelled surface electrode for direct connection

Order #	Size	Lead Length	Instrument End
9013L0453	9 x 6 mm	50 cm	1.5 mm female TP

Each package contains 12 electrodes

Tab electrode - to be used with connection cable

Order #	Size	Lead Length	Instrument End
9013S0232	10 x 10 mm	n/a	n/a

Each package contains 100 electrodes

Connection cable for Tab electrodes

Order #	Color	Lead Length	Instrument End
9013P0602	Green	80 cm	1.5 mm female TP
9013P0642	Red	80 cm	1.5 mm female TP

Each package contains 2 cables



Benefits of a large system in a "smartphone-size" device

Clavis provides muscle and nerve localization wherever needed to accurately guide and monitor injection of medicinal products. It provides EMG recording with audio and stimulation ranging from 0-15 mA with speed and reliability while efficiently operating on battery power.

Guidance for injection of medicinal products where & when you need it

Clavis is quick and easy to operate. Provides guidance for accurate injection of medicinal products at the exact point of treatment.

This unique device takes up minimal space, so medical professionals can get in and out of small spaces quickly.

Instantly switch between stimulation and recording

Single-handedly switch between stimulation and recording functions.

When seconds really count, Clavis' fast and simple operation helps to save time.




Dantec Clavis & Bo-ject needles – made for each other

Bo-ject needles are:

- Designed to minimize patient discomfort by applying a thin trocar-shaped construction and low friction coating to ease insertion
- Compatible with all types of reference electrodes and fits all types of syringes with Luer-lok™ or slip fit connectors
- Easy to attach to Clavis and ready to use quickly
- Ergonomically designed and lightweight for greater control

To ensure consistently reliable results, Clavis is designed to work with Bo-ject hypodermic needle electrodes.

These needles have been specifically developed and cleared for injection of medicinal products such as botulinum toxin.



Service

Natus Neurology is committed to providing exemplary service to our customers. Our dedicated and experienced Customer Service Team will assist with every aspect of an order. To support our products, we provide factory-trained Field Technicians and Clinical Application Specialists for onsite support. Additionally, we provide an in-house Technical Support Team, staffed with experts, and a strong distribution network in International Markets to offer a wide range of service options. Allowing our customers more time to care for their patients is our goal. Customer loyalty is our reward.

Supplies

Natus Neurology offers a full range of neurodiagnostic accessories and supplies promoting patient comfort. Our dedicated customer service team provides a streamlined order and shipping process to save you time and money.

To learn more about Natus Neurology Service Programs or our full line of Supplies and Accessories, contact your local distributor or sales representative.
US Customers Call: **1-800-356-0007**



Natus Neurology Incorporated

3150 Pleasant View Road
Middleton, WI 53562 USA
Tel: 1-800-356-0007
1-608-829-8500
Fax: 1-608-829-8709
www.natus.com

Natus Europe GmbH

Robert-Koch-Str. 1
82152 Planegg
Germany
Tel: +49 (0)89 / 83942-0
Fax: +49 (0)89 / 83942-777

© 2014 Natus Medical Incorporated. All Rights Reserved. All product names appearing on this document are trademarks or registered trademarks owned, licensed to, promoted or distributed by Natus Medical Incorporated, its subsidiaries or affiliates.